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Make Research Easier

Description

This kit uses a negative selection method to isolate CD3 T cells from human PBMCs through column-free magnet. During the separation process, the antibodies and Nanobeads do not contact with the CD3 T cells, allowing for the preservation of the T cells in their most original state.

- Easy to operate
- Purity up to 99%
- Fast, only 18-20 minutes for separation

The kit labels non-CD3 T cells using biotin-conjugated antibodies, followed by incubation with streptavidin-coated magnetic nanoparticles for coupling. The separation is then performed using a magnetic field without column. Desired cells can simply be transferred to a new tube for separation, and the isolated cells can be immediately used for subsequent applications, such as flow cytometry, cell culture, and functional assays.

Component Descriptions

COMPONENT NAME	STORAGE	BUFFER	FUNCTION
EasyIso™ Human CD3 T cell Isolation Cocktail	Store at 4°C. Do not freeze.	PBS, 0.1% BSA	The Isolation Cocktail identifies non-CD3 T cells, marking the non-target cells.
EasyIso™ Streptavidin Nanobeads	Store at 4°C. Do not freeze.	PBS, 0.1% BSA	Streptavidin-coated magnetic nanoparticles bind to non-CD3 T cells, which are then removed using a magnetic rack.

BSA - bovine serum albumin; PBS - phosphate-buffered saline

Components should be shipped at 4°C.

Preparation

- Buffer: FACS Buffer (1 x PBS+0.02% EDTA+2% FBS)
- 70 µm sterile nylon filter; sterile centrifuge tubes
- Centrifuge, Magnetic rack

Note

- It is recommended to use lymphocyte separation medium to isolate PBMCs; red blood cell lysis buffer is not recommended.
- Use FACS buffer instead of PBS during operations to improve cells yield.
- Ensure the osmotic pressure of the FACS buffer is stable, as immune cells are sensitive.
- Maintain sterile techniques at every step to prevent contamination.

Operating Steps

STEPS	INSTRUCTIONS	DOSAGE and TIME
1	Prepare the PBMC single-cell suspension at the indicated concentration	1×10^8 cells/mL
2	Add EasyIso™ Human CD3 T cell Isolation Cocktail to the sample	100 μ L/mL
3	Mix and incubate	RT for 10 min
4	Vortex the Streptavidin Nanobeads	Vortex for 30 s
5	Add Streptavidin Nanobeads to the sample	100 μ L/mL
6	Mix and incubate	RT for 5 min
7	Add FACS buffer to the sample	Top up to 2.5 mL (5 mL tube) Top up to 7.5 mL (15 mL tube)
8	Place the tube into the magnet and incubate	RT for 5 min
9	Pick up the magnet, transfer the sample into a new collection tube to collect the sample	Isolated cells are ready for use

Data


Figure. Using the Immuno-X™ series products, CD3 T cells are isolated and activated. (A) The CD3 T cells is isolated from healthy human PBMCs, achieving a purity of up to 98.7% with CD3 flow antibodies for staining (#SN3101). (B) CD3 T cells are labeled with CFSE, followed by stimulation using Human ActBeads™ T cell Activator CD3/CD28 (#AH2001). After 3 days, flow cytometric analysis reveals proliferation of T cells

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Eynos LifeScience co., Ltd